

# SkinPen<sup>®</sup>

by Bellus Medical



A microneedling device indicated to treat melasma,  
acne scarring and surgical scars.

AN ARTG INCLUDED MEDICAL DEVICE  
NOW IN AUSTRALIA AND NEW ZEALAND

## The SkinPen Experience

As we age the natural architecture of our skin deteriorates, which leads to fine lines, wrinkles, large pores and uneven skin texture.

Microneedling is a way to address these signs of aging.

It creates thousands of micro-injuries in the skin in an effort to trigger the body's natural wound healing process.

This assists in the promotion of collagen and elastin, the two proteins predominantly responsible for the structure of the skin.

## The SkinPen Results



### FEMALE, AGE 22 | 6 SKINPEN PROCEDURES\*

Photos courtesy of: Cathy A. Presnick, Licensed Aesthetician; A Perfect Complexion, LLC. Melbourne, FL.



### FEMALE, AGE 51 | 2 SKINPEN PROCEDURES\*

Photos courtesy of: Heather Dawson, NP-C; Dawson Wellness, McKinney, TX



### FEMALE, AGE 38 | 5 SKINPEN PROCEDURES\*

Photos courtesy of: Dr. Christie Matter; North Texas Dermatology, Plano, TX



### FEMALE, AGE 45 | 3 SKINPEN PROCEDURES\*

Photos courtesy of: Elizabeth Steigner, Medical Aesthetician; Turpin, MD Plastic Surgery & Medical Aesthetics, Orange, CA

## AN ARTG INCLUDED MEDICAL DEVICE #297614 AND #296066

Years of Bellus Medical's proactive work in design, development and studies, has resulted in a Class IIa medical device included on the ARTG. SkinPen effectively treats melasma, acne scarring and surgical scar.

### DESIGNED AND ENGINEERED FOR EXCELLENCE IN EVERY STEP

**INDUCTIVE CHARGING** - Inductive battery charging allows for 100% charge to be achieved between every patient.

**CORDLESS** - Cordless technology improves mobility.

**POWER SWITCH WITH AUTO-SHUT OFF** - The device is designed to automatically switch off if the correct needle length cannot be achieved; tolerance is up to 0.001mm

**LED SMART TECHNOLOGY** - Battery life indicator with LED technology

**BATTERY** - ROHS compliant

**MOTOR WITH INTUITIVE SPEED** - Intuitive Speed Control regardless of skin thickness or condition, motor output is automatically adjusted.

**ERGONOMIC DESIGN** - The SkinPen® ergonomic design increases treatment comfort for ease of use

**BIOSLEEVE** - Protects SkinPen during treatment

**ULTRASONIC SEAL** - Ultrasonic seal between pen and cartridge

**PRECISION ADVANCED CARTRIDGE UNIT™** - The Precision Advanced Cartridge Unit (ACU) is a single use cartridge ensuring long term reliability of the pen  
ARTG 297614

**14 STAINLESS STEEL NEEDLES, SMOOTH & STERILE** - The Precision Advanced Cartridge Unit contains 14 smooth, stainless steel, needles (~250 um diameter)

**PRECISION DEPTH DIAL**

**INTERNAL RECIPROCATING MECHANISM**

**SINGLE USE LOCKOUT FEATURE** - The Precision Advanced Cartridge Unit (ACU) is sterile with a Lockout feature for single use only, ensuring the long term accuracy of the pen.

ARTG 296066

## Patient Treatment Kit



### KIT INCLUDES:

SkinPen Precision Advanced Cartridge Unit (ACU), Skinfuse Lift HG, Skinfuse Rescue Calming Complex & Disposable Bio-Sheath

ARTG 296065

Exclusively distributed in Australia and New Zealand by



\*Cases typical of results at Bellus Medical, individual results may vary

Exclusively distributed in Australia  
and New Zealand by



---

For more information please contact your AST Account Manager or call  
**AU** [customerservice@advskin.com.au](mailto:customerservice@advskin.com.au) | **NZ** [customerservice@advskin.co.nz](mailto:customerservice@advskin.co.nz)  
**AU** 1800 648 851 | **NZ** 0800 238 754

---

**SkinPen**<sup>®</sup>  
by Bellus Medical

**BELLUS**  
**MEDICAL**